



**AUTHORIZED
DISTRIBUTOR**

2026

**PRODUCT
CATALOG**





NEURO PATTIES

RAY-COT®	2
DELICOT®	4
TAPERED DELICOT®	5
UNIQCOT™	6
TELFA®	7
AMERICOT®	8

SINGLE-USE DEVICES

BIPOLAR FORCEPS	10
RHINOSTOP™ NASAL PACKING	12


OTHER SPECIALTIES

COTTON BALLS	13
OPHTHALMIC PLEDGETS	13

RAY-COT[®] Rayon Patties

AN IMPROVED VERSION OF THE STANDARD RAYON NEUROSURGICAL PATTIE WITH SUPERIOR STRENGTH, ABSORBENCY, AND PERFORMANCE.

- Single x-ray detectable barium line allows for superior pliability and softness
- Smooth precise edges reducing the potential for loose fibers
- Thin locator string that has been sewn in for increased safety
- Packaged in quantities of 200 per box (20 pouches of 10)



#60-02

6 mm x 6 mm
¼ in x ¼ in



#60-04

6 mm x 76 mm
¼ in x 3 in



#60-05

10 mm x 76 mm
⅜ in x 3 in



#60-06

13 mm x 13 mm
½ in x ½ in



#60-07

13 mm x 25 mm
½ in x 1 in



#60-08

13 mm x 38 mm
½ in x 1½ in



#60-09

13 mm x 76 mm
½ in x 3 in



#60-10

19 mm x 19 mm
¾ in x ¾ in




#60-12

25 mm x 25 mm
1 in x 1 in



#60-13

25 mm x 76 mm
1 in x 3 in



#60-14

38 mm x 38 mm
1½ in x 1½ in



#60-15

76 mm x 76 mm
3 in x 3 in



#60-16

6 mm x 38 mm
¼ in x 1½ in



#60-22

38 mm x 76 mm
1½ in x 3 in





#60-29

13 mm x 51 mm
½ in x 2 in

RAY-COT® Rayon Strips

AN IMPROVED VERSION OF THE STANDARD RAYON NEUROSURGICAL PATTIE WITH SUPERIOR STRENGTH, ABSORBENCY, AND PERFORMANCE


- Single x-ray detectable barium line allows for superior pliability and softness
- Strips come without a locator string
- Packaged in quantities of 200 per box (20 pouches of 10)



DELICOT[®] Rayon Patties

ULTRA-THIN NEUROSURGICAL PATTIES DESIGNED FOR SKULL BASE KEYHOLE MICROSURGERY, FEATURING A 0.33 MM PROFILE THAT OFFERS EXCEPTIONAL PLIABILITY AND PRECISE MOLDING TO WRAP DELICATE STRUCTURES, PROTECT TISSUE, AND MAINTAIN CLEAR VISIBILITY UNDER THE MICROSCOPE.

- Thinnest neurosurgical pattie available
- Single x-ray detectable barium line allows for superior pliability and softness
- Smooth precise edges reducing the potential for loose fibers
- Thin locator string that has been sewn in for increased safety
- Packaged in quantities of 200 per box (20 pouches of 10)*



#63x00

2 mm x 20 mm
3/16 in x 25/32 in



#63-01

5 mm x 8 mm
3/16 in x 5/16 in



#63-02

5 mm x 13 mm
3/16 in x 1/2 in



#63-03

10 mm x 13 mm
3/8 in x 1/2 in



#63-04

15 mm x 20 mm
9/16 in x 13/16 in



#63-05

4 mm x 50 mm
1/2 in x 2 in



#63-06

6 mm x 50 mm
1/4 in x 2 in



#63-07

9 mm x 50 mm
11/32 in x 2 in



#63-08

13 mm x 60 mm
1/2 in x 2 3/8 in



#63-09

20 mm x 60 mm
13/16 in x 2 3/8 in



#63-10

30 mm x 70 mm
1 3/16 in x 2 3/4 in

*63x00 Packaged in quantities of 100 per box (10 pouches of 10)


DELICOT® Tapered Rayon Patties

A COMBINATION OF UNPARALLELED FLEXIBILITY, THINNESS, AND TAPERED SHAPE OF THIS CLASS OF PATTIES PROVIDES THE SURGEON WITH GREAT ADVANTAGES WORKING WITHIN SMALL OPERATIVE SPACES SEAMLESSLY


- Thinnest neurosurgical pattie available
- Single x-ray detectable barium line allows for superior pliability and softness
- Smooth precise edges reducing the potential for loose fibers
- Thin locator string that has been sewn in for increased safety
- Packaged in quantities of 200 per box (20 pouches of 10)




#63-07T
4 mm x 9 mm x 48 mm



#63-09T
8 mm x 19 mm x 64 mm




#63-12T
8 mm x 19 mm x 76 mm




#63-08T
6 mm x 19 mm x 51 mm



#63-13T
6 mm x 13 mm x 9 mm



#63-14T
19 mm x 25 mm x 25 mm




#63-15T
19 mm x 25 mm x 76 mm

UNIQCOT™ Rayon Patties & Strips

THE GOLDEN MEAN BETWEEN RAY-COT® AND DELICOT®, FOR THOSE LOOKING FOR A THINNER PATTIE FOR BETTER VISIBILITY WITHOUT COMPROMISING ABSORBENCY.

- Same trusted properties of Ray-Cot but only 0.5 mm thick
- Single x-ray detectable barium line allows for superior pliability and softness
- Smooth precise edges reducing the potential for loose fibers
- Thin locator string that has been sewn in for increased safety
- Packaged in quantities of 200 per box (20 pouches of 10)



#67-02
6 mm x 6 mm
¼ in x ¼ in



#67-04
6 mm x 76 mm
¼ in x 3 in



#67-06
13 mm x 13 mm
½ in x ½ in



#67-07
13 mm x 25 mm
½ in x 1 in



#67-08
13 mm x 38 mm
½ in x 1½ in



#67-09
13 mm x 76 mm
½ in x 3 in



#67-10
19 mm x 19 mm
¾ in x ¾ in



#67-12
25 mm x 25 mm
1 in x 1 in



#67-13
25 mm x 76 mm
1 in x 3 in



#67-29
13 mm x 51 mm
½ in x 2 in

UNIQCOT™ Strips
No locator string



#67-81
6 mm x 152 mm
¼ in x 6 in



#67-82
13 mm x 152 mm
½ in x 6 in



#67-84
25 mm x 152 mm
1 in x 6 in

TELFA[®] Surgical Patties & Strips

PERFORATED NATURAL POLYMER LAYERS SURROUND A COTTON CORE, ALLOWING FLUID TO PASS THROUGH WHILE MINIMIZING SURFACE ADHESION AND MAINTAINING GENTLE TISSUE CONTACT

- Telfa[®] outer layer surrounding an absorbent cotton layer
- Remarkably low tissue adhesion with high absorbency
- Single x-ray detectable barium line allows for superior pliability and softness
- Thin locator string that has been sewn in for increased safety
- Packaged in quantities of 200 per box (20 pouches of 10)



Telfa[®] is a registered trademark of KPR U.S., LLC



#80-04
6 mm x 76 mm
¼ in x 3 in



#80-05
10 mm x 76 mm
¾ in x 3 in



#80-06
13 mm x 13 mm
½ in x ½ in




#80-08
13 mm x 38 mm
½ in x 1½ in




#80-09
13 mm x 76 mm
½ in x 3 in



#80-11
19 mm x 76 mm
¾ in x 3 in



#80-12
25 mm x 25 mm
1 in x 1 in



#80-14
76 mm x 76 mm
3 in x 3 in




#80-13
25 mm x 76 mm
1 in x 3 in

TELFA[®] Strips
No locator string



#81-06
13 mm x 152 mm
½ in x 6 in




#81-13
25 mm x 152 mm
1 in x 6 in

AMERICOT[®] Cotton Patties

HIGHLY PLIABLE AND SOFT COTTON PATTIE THAT ALLOWS FOR UP TO 11 TIMES ITS DRY WEIGHT TO BE ABSORBED WHEN FULLY SATURATED

- Precision cut providing smooth edges, reducing the potential for loose fibers
- Single x-ray detectable barium line allows for superior pliability and softness
- Thin locator string that has been sewn in for increased safety
- Packaged in quantities of 200 per box (20 pouches of 10)



#20-06s

13 mm x 38 mm
½ in x 1½ in



#20-05s

13 mm x 25 mm
½ in x 1 in



#20-04s

13 mm x 13 mm
½ in x ½ in



#20-01s

6 mm x 6 mm
¼ in x ¼ in



#20-19s

3 mm x 3 mm
¼ in x ¼ in



#20-07s

13 mm x 76 mm
½ in x 3 in



#20-22s

10 mm x 10 mm
¾ in x ¾ in



#20-08s

19 mm x 19 mm
¾ in x ¾ in



#20-09s

25 mm x 25 mm
1 in x 1 in



#20-10s

25 mm x 76 mm
1 in x 3 in



#20-13s

19 mm x 76 mm
¾ in x 3 in



#20-12s

76 mm x 76 mm
3 in x 3 in



#20-14s

10 mm x 19 mm
¾ in x ¾ in



#20-27s

19 mm x 51 mm
¾ in x 2 in



#20-25s

10 mm x 25 mm
¾ in x 1 in



#20-26s

10 mm x 51 mm
¾ in x 2 in



#20-32s

13 mm x 51 mm
½ in x 2 in




#20-30s

6 mm x 38 mm
¼ in x 1½ in

AMERICOT[®] Cotton Strips

AMERICOT[®] IS CAPABLE OF HOLDING GREATER THAN 11 TIMES ITS DRY WEIGHT IN FLUID WHEN FULLY SATURATED

- Precision cut providing smooth edges, reducing the potential for loose fibers
- High degree of wet strength and conformity
- Single x-ray detectable barium line allows for superior pliability and softness
- Strips come without a locator string
- Packaged in quantities of 200 per box (20 pouches of 10)



#21-01

6 mm x 152 mm
¼ in x 6 in



#21-02

13 mm x 152 mm
½ in x 6 in




#21-04

25 mm x 152 mm
1 in x 6 in

ATLAS™ CHOICE Ultra Non-stick Bipolar Forceps

THE ATLAS™ CHOICE DISPOSABLE NON-STICK BIPOLAR FORCEPS FEATURE TIPS CREATED SPECIFICALLY FOR MICROSURGERY (MICRO TIPS) AND TUMOR AND SPINE (MACRO TIPS) PROCEDURES.

- Solid silver precision tips
- Stainless steel body for improved durability
- Low-profile slim tines
- Dual irrigation options
- Packaged in quantities of 1 per box



Length

Tip Size and Product Codes

8"

0.5 mm Micro Tip

Non-Irrigating: BIP-ATMI-04

Irrigating: BIP-ATMI-04-IR

1.0 mm Micro Tip

Non-Irrigating: BIP-ATMI-05

Irrigating: BIP-ATMI-05-IR

1.5 mm Macro Tip

Non-Irrigating: BIP-ATMA-03

Irrigating: BIP-ATMA-03-IR

2.0 mm Macro Tip

Non-Irrigating: BIP-ATMA-04

Irrigating: BIP-ATMA-04-IR

9"

0.5 mm Micro Tip

Non-Irrigating: BIP-ATMI-07

Irrigating: BIP-ATMI-07-IR

1.0 mm Micro Tip


Non-Irrigating: BIP-ATMI-08

Irrigating: BIP-ATMI-08-IR

1.5 mm Macro Tip


Non-Irrigating: BIP-ATMI-09

Irrigating: BIP-ATMI-09-IR



Superior Tip Functionality

Our solid silver tips create a superior non-stick effect. Perfect tip alignment provides balance and control for seamless dissection without a scissoring effect.



Advanced Heat Transport

Advanced heat transport technology efficiently regulates tine and tip temperature to minimize any tissue damage or charring.

RhinoStop™ Nasal Packing

A NASAL PACKING SOLUTION FOR EFFECTIVE EPISTAXIS MANAGEMENT. MADE FROM ULTRA-CLEAN POLYVINYL ALCOHOL (PVA) USING PROPRIETARY AIR-FOAMING TECHNOLOGY TO REDUCE MICROBIAL GROWTH AND CHEMICAL RESIDUE. ITS SOFT, SUPPORTIVE STRUCTURE ABSORBS QUICKLY, MINIMIZES ADHESION, AND REDUCES IRRITATION FOR COMFORTABLE, RELIABLE HEMOSTASIS.

- Fast absorption of blood and fluids
- Lint- and fiber-free material
- Ultra-clean process removes chemical residue
- Soft and supportive for patient comfort
- Proprietary air-foaming technology

Standard Type



Sinus Type



Epistaxis Type



Standard Type

SKU	SIZE	STRING/AIRWAY TUBE
140-01	45x20x10	Yes/No
140-02	45x20x15	No/No
140-03	45x20x15	Yes/No
140-04	45x20x15	Yes/Yes
140-05	80x20x10	No/No
140-06	80x20x10	Yes/No
140-07	80x20x15	No/No
140-08	80x20x15	Yes/No
140-13	80x20x15	Yes/Yes

Epistaxis Type


SKU	SIZE	STRING
140-09	100x25x15	No
140-10	100x25x15	Yes

Sinus Type


SKU	SIZE	STRING
140-11	35x12x6	Yes
140-12	35x12x9	Yes

X-RAY DETECTABLE Cotton Balls


- 100% American Cotton
- Each premium strung cotton ball contains an x-ray detectable marker for increased safety and accountability
- Double strung locator string allow for easier traceability post-procedure
- Packaged in quantities of 100 per box (20 pouches of 5)



X-Large
#30-05
32 mm
1¼ in
(Strung)



Large
#30-04
25 mm
1 in
(Strung)



Medium
#30-03
19 mm
¾ in
(Strung)



Small
#30-00
6 mm
¼ in
(Strung)



X-Small
#30-07
3 mm
⅛ in
(Strung)



X-Small
#31-07
3 mm
⅛ in
(Non-Strung)



Small
#31-00
6 mm
¼ in
(Non-Strung)

COMFORTPACK™ Ophthalmic Pledgets

THE COMFORTPACK™ EYE PLEDGETS FOR OPHTHALMOLOGY

- Absorbs anesthetic solutions for precise and effective numbing
- Eliminates waste and increases cost savings
- Sterile double strung cotton pledgets are a multi-use product for the clinic and operating room

#110-01
ComfortPack™ Kit
> Two 5/32" (4 mm) double strung pledgets
> One dispensing cup and swab
> Packaged in quantities of 60 (30 pouches of 2)

#111-01
Ophthalmic Pledgets 5/32"
> Five 5/32" (4 mm) double strung pledgets only per pouch
> Packaged in quantities of 100 (20 pouches of 5)

#111-02
Ophthalmic Pledgets 7/32"
> Five 7/32" (5.5 mm) double strung pledgets only per pouch
> Packaged in quantities of 100 (20 pouches of 5)



**Sterile Mate Is Now On
The AskIsaac.io
intelligence
platform....the future of
medical supply chain
fulfillment and facilities
Management Solutions**

Order Now

1-800-988-9832 • 1-281-456-3458

Orders@sterilemate.com

www.sterilemate.com

